



SoGo[®] Medicinal Stain Removal System

The SoGo Medicinal Stain Removal System is a linen rewash system designed to return linens stained with Betadyne or iodine to working inventory without damaging the linen tensile strength or color. The SoGo Medicinal Stain Removal System contains SoGo "Oil Break" and SoGo "S.M.R."

SoGo "Oil Break" is a liquid laundry detergent / stain removal agent, specially designed to be used in a laundry cycle with SoGo "S.M.R.". SoGo "Oil Break" works synergistically with SoGo "S.M.R." to attack medicinal stains.

SoGo "S.M.R." is a powdered laundry stain remover whose unique, controlled bleaching action safely removes residual medicinal stains as well as mildew.

DIRECTIONS

Load machine to capacity, but do not overload.

STEP 1	185°F ¹ Water at High Level	Flush	Run for 3 minutes
STEP 2	185°F* Water at Low Level	Add 16 oz. SoGo "Oil Break" plus 32 oz. SoGo "S.M.R." for each 100 lbs. of rated machine capacity.	Run for 30 minutes
STEP 3	Follow with four hot, high level rinses, then sour to finish. Intermediate extracts between the rinses are recommended.		

¹ Hotter water will yield better stain removal results. At temperatures below 165°F, increase the running time in step two to 45 minutes. Minimum temperature is 140°F.

NOTE: SoGo Medicinal Stain Removal System can be used on whites and color-fast colors. Add products manually after temperature level is reached. Quantities indicated are for horizontal axis commercial laundry machines only.

PRODUCT OVERVIEW